




Dip the **DrugControl PHENCYCLIDINE** strip into urine specimen



Pipette 3 – 4 drops of urine specimen onto the adsorbent tip of **DrugControl PHENCYCLIDINE** strip

DIRECTIONS FOR USE

- 1 Allow the test, urine specimen, and / or controls to reach room temperature (15 – 30 °C) prior to testing.
- 2 Bring the pouch to room temperature before opening it.
- 3 Remove the test device from the sealed pouch and use it as soon as possible.
- 4 Place the test device on a clean and level surface.
- 5 Add 3 – 4 drops of urine (approx. 240µl) to the specimen well or immerse the test vertically in the urine specimen for at least 10-15 seconds.
- 6 Do not pass the maximum line (Dipping line) on the test when immersing the strip.
- 7 Place the test on a non-absorbent flat surface, start the timer and wait for the red line(s) to appear.
- 8 The result should be read at 10 minutes.

INTERPRETATION OF RESULTS



Negative



Positive



Invalid



Invalid

- Negative:*** Two lines appear. One red line should be in the control region (C), and another apparent red or pink line should be in the test region (T). This negative result indicates that the Phencyclidine concentration is below the detectable level (25 ng/mL).
- Positive:** One red line appears in the control region (C). No line appears in the test region (T). This positive result indicates that the Phencyclidine concentration exceeds the detectable level (25 ng/mL).
- Invalid:** Control line fails to appear. Insufficient specimen volume or incorrect procedural techniques are the most likely reasons for control line failure. Review the procedure and repeat the test using a new test device. If the problem persists, discontinue using the lot immediately and contact distributor / manufacturer.

* **Note:** The shade of red in the test line region (T) may vary, but it should be considered negative whenever there is even a faint pink line.

Detailed informations regarding this lateral flow chromatographic immunoassay are given in the respective instructions for use:

PHENCYCLIDINE TEST

